Multi-Safe® Cu 375 short

Intrauterinapparat
Intrauterine Device
Dispositif Intra-Utérin
Dispositivo Intrauterino

Multi-Safe® Cu 375 short

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GB

Multi-Safe® Cu 375 short

Intrauterine Contraceptive Device

PHYSICIANS'S INSTRUCTIONS

Presentation

The Multi-Safe Cu 375 short Intrauterine Device (IUD) is made of polyethylene. The vertical arm is wound with copper wire. The copper surface is 375 mm². Dimension: 19,5 mm ± 1 mm wide and 29,4 mm ± 0,2 mm long. The wire is radiopaque.

Characteristics

Multi-Safe Cu 375 short IUD offers almost complete protection against pregnancy. It remains effective longer than any other non-copper IUD-based reversible family planning method. Correctly inserted, the Multi-Safe Cu 375 short IUD is safe for women at low risk of sexually transmitted disease. The Multi-Safe Cu 375 short is an intrauterine contraceptive device made of polyethylene. The side arms are flexible and shaped in such a way as to keep the IUD adjacent to the fundus, without stretching the uterine cavity or touching the entrance of the fallopian tubes.

Multi-Safe Cu 375 short IUD prevents pregnancy by blocking fertilization.

The contraceptive effect is likely due to a pronounced sterile inflammatory reaction, which takes place as a result of a foreign body response in the uterus. The concentration of various types of white blood cells, prostaglandins and enzymes in uterine and tubal fluids increases markedly. These changes interfere with the movement of sperm in the genital tract, reducing their potential fertilizing capacity, so that fertilization is not possible. In the fallopian tubes, where fertilization is thought to take place, there are found fewer sperm in copper IUD users than in non-users. Thus, the primary action is most likely altering the function or survival of sperm and ova before they can meet.

Indication

Female contraception (long lasting)

Contraindications (absolute)

- Malignant diseases of the genital tract
- Vaginal bleeding
- Pregnancy
- Past history of ectopic pregnancy of predisposing factors
- Infections of the genital tract
- Sexually transmitted diseases during the last 12 months (except bacterial vaginitis, repeated herpes infection, hepatitis B)
- Abortion with infection during the last 3 months, pelvic inflammatory disease (PID)
- Uterine malformations (congenital or acquired)
- Allergy to copper

Contraindications (relative)

- Valvular heart disease
- Anaemia
- Coagulation disorders
- Anti-inflammatory treatment
- Wilson's disease
- Multiple exposures to different sexual partners
- Nulliparity

Precautions

Before inserting the IUD, a thorough medical history and an examination of the pelvic and abdominal cavity as well as a cervical smear are mandatory. Pregnancy, genital infections or sexually transmitted diseases have to be excluded. The position of the uterus and the size of the uterine cavity must be determined to ensure correct insertion of the IUD. Multi-Safe Cu 375 short is designed for women with an uterine cavity depth of 6 - 9 cm. One month after the insertion of an IUD the woman must be re-examined to determine whether the IUD is properly placed and if there are signs of infection. Subsequent follow-up examinations are to be performed annually, if clinically indicated even more frequently. Pelvic inflammatory disease during IUD use should be treated without delay. For this reason the user must be instructed to report to

her physician if there are suspicious signs. This can be ascertained by gynaecological examination and/or ultrasound (if available). If a treatment with antibiotics is not effective after 48 hours or if there are signs of PID, the IUD must be removed immediately.

Excessive bleeding or dysmenorrhoea during the first cycle after insertion should also be carefully examined to ascertain if this is caused by the IUD, in which case it might have to be removed. The possibility of perforation of the uterus during the insertion should always be considered, especially if the nylon threads are invisible or cannot be pulled out of the cervical canal. If there are any doubts about the position of the IUD (e.g. if the IUD is extremely difficult or painful to insert) the appropriate diagnostic techniques should be used (flat x-ray of the pelvis, ultrasonography, hysteroscopy, laparoscopy).

If the nylon threads appear to be longer than at insertion, an ultrasonogram should be carried out to determine if the IUD has been displaced, what might decrease its contraceptive efficacy. If the threads aren't visible in front of the cervix at a follow-up examination, at first pregnancy must be excluded. The threads may have possibly been drawn up into the cervical canal or the uterus and usually reappear during the next menstruation. In the case of doubt the IUD can be identified by feeling carefully with a suitable instrument, by ultrasound diagnostics or by an x-ray after exclusion of a pregnancy.

If IUD users get accidentally pregnant, it is mandatory to determine (by ultrasound) whether this is an intrauterine or ectopic pregnancy. Up to the end of the first trimester, if the threads are visible, the IUD should be removed. After that the patient should be offered the option of elective abortion as soon as possible, since there is an increased risk of PID and many other obstetrical problems, e.g. premature labour, placenta praevia and ablatio placentae. However, if the woman wishes to continue her pregnancy, close monitoring is mandatory. She should be informed about the risks of keeping the IUD in situ.

Perforation

Perforation or penetration of the uterine corpus or cervix by the IUD may occur, most often during insertion. The risk is increased during the postpartum period and slightly increased if the woman is fully breast feeding. Such a device must be removed as soon as possible.

Interactions

The available experience indicates that, in general, drugs interfering with the contraceptive efficacy of Multi-Safe Cu 375 short are highly unlikely. However, published reports appear to show diminished efficacy in the presence of long term use of non-steroidal anti-inflammatory drugs (especially acetyl salicylic acid) and of corticoids. Short term use in the treatment of dysmenorrhoea with non-steroidal anti-inflammatory drugs does not appear to reduce contraceptive efficacy.

Side effects

During use of IUDs stronger menstrual bleeding, slight intermenstrual bleeding, anaemia, dysmenorrhoea, pain in the lower abdomen or sacral area as well as lower abdominal infections may occur. The risk of an ectopic pregnancy in the event of failure of the method is increased. The IUD as well as parts of the same may penetrate into the wall of the uterus or perforate it. Allergic skin reactions may occur.

Physician's instructions to IUD users

Users should be informed in detail as to the advantages and disadvantages of IUD contraception, not only so that they understand how it works, but above all so that any complications can be detected early. The user must learn how to feel the threads emerging from the cervical canal.

Insertion

First eliminate the possibility of pregnancy!

The best time for insertion is during menstruation to prevent insertion during non-diagnosed pregnancy. At this time the external and internal cervical os are physiologically dilated. This facilitates the insertion of the IUD without the need to dilate the canal in most instances. Alternatively, the IUD may be introduced within 5 days of unprotected coitus or within 15 minutes of delivery of the placenta or abortion in the first trimester. In these last two cases a higher expulsion rate must be accepted. Insertion immediately after unprotected coitus can increase the risk of PID.

If the IUD cannot be inserted immediately after delivery of the placenta or abortion, insertion should be delayed for at least six weeks.

In case of caesarean section insertion should be delayed for 12 weeks after delivery, although this is controversial.

Prior to insertion the vagina, cervix and cervical canal should be cleansed with an antiseptic solution, using e.g. a wooden cotton bud. It is essential to determine the exact position of the uterus by bimanual pelvic palpation so that the Multi-Safe Cu 375 short can be inserted along its longitudinal axis. This can be accomplished by grasping the anterior or posterior lip of the cervix, depending on whether the uterus is anteverted or retroverted.

In case of vasovagal reactions after the use of a forceps a local anaesthetic can be injected in and around the cervix. Hysterometry should be carried out and the ascertained depth of the uterus be marked on the Multi-Safe Cu 375 short inserter tube with the aid of the movable flange.

Step-by-Step loading instructions

The IUD may be used by trained medical staff only. In order to minimize the risk of contamination we recommend using sterile gloves.

- Open the sterile packaging of the Multi-Safe Cu 375 short half-way and mark the previously sounded length of the uterus with the movable flange around the insertion tube.


fig.1

- Remove the loaded insertion tube from the package and insert it gently into the uterus until the movable flange touches the cervical os which means that the tip of the IUD has now reached the fundus.


fig.2

- Pull off the insertion tube releasing the Multi-Safe Cu 375 short in the uterus.


fig.3

Removal

Multi-Safe Cu 375 short should be replaced after five years.

Again, during menstruation is the most appropriate time to remove the IUD, since both the internal and external cervical os are fully dilated. Grasp the threads of the IUD with a tenaculum and pull it along the longitudinal axis of the uterus. Try to insert the forceps at the entrance to the cervical canal in order to grasp the vertical arm of the Multi-Safe Cu 375 short as soon as it has passed the internal cervical os. This avoids excessive tension on the threads, which could cause them to break. While doing this, distract the patient's attention by asking her to cough, and then remove the Multi-Safe Cu 375 short with one firm tug.

Insertion and removal of the IUD could cause slight pain and bleeding. The intervention may also precipitate a vasovagal attack or an epileptic seizure.

After removal is complete, inspect the Multi-Safe Cu 375 short to see that none of its arms have been left in the uterine cavity. In case of a lost IUD in the uterine cavity either hysteroscopy or ultrasonography and x-rays should be used to determine its location; curettage may be advisable. In very rare cases of uterine perforation laparoscopy may be needed.

Packaging

1 x 1 sterile Multi-Safe Cu 375 short
10 x 1 sterile Multi-Safe Cu 375 short

Pharmaceutical information

Each IUD is sterilized with ethylene oxide and is intended for single use only. Do not resterilize. Do not use if package is damaged or open. Do not use after the expiry date printed on the sterile package. After removal the Multi-Safe Cu 375 short is to be disposed of in accordance with the guidelines in force for dealing with potentially infectious material.

Shelf-life: Four years.

Storage: Store in a dry place, below 35 °C. Protect from direct sunlight and moisture.

List of excipients: Copper, Polyethylene, Polyamide 6.

Incompatibilities: No known incompatibilities.

Nature and contents of container: The device with accessories has been packed in a heat sealed sterilised pouch made of Tyvek + PET/PE.

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